Discotom-6

Instruction Manual

Manual No.: 15677001

Date of Release 03.02.2006

Discotom-6
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	15
Quick Reference Guide.....	27

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2006.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Discotom-6

Safety Precaution Sheet

To be read carefully before use

1. The operator(s) should be fully instructed in the use of the machine and its cut-off wheels according to the Instruction Manual and the instructions for the cut-off wheels.
2. The machine must be placed on a safe and stable support table. All safety functions and guards of the machine must be in working order.
3. Use only intact cut-off wheels. The cut-off wheels must be approved for min. 50 m/s.
4. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
5. The workpiece must be securely fixed in the quick-clamping device or similar. Large or sharp workpieces must be handled in a safe way.
6. The cutting handle should be lowered slowly and carefully, in order to avoid breaking the cut-off wheel.
7. Struers recommends the use of an exhaust system as the materials being cut may emit harmful gasses or dust.
8. The rubber curtains on the left hand side make it possible to cut long workpieces with a height of up to 80mm.
9. Do not put your hand through the rubber curtain during operation of the machine.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.)

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents	2
Recirculation Cooling Unit.....	2
Unpacking Discotom-6.....	2
Placing Discotom-6	2
Getting Acquainted with Discotom-6.....	3
Discotom-6 Rear View	4
Noise Level.....	4
Power Supply	4
Direction of the Cut-off Wheel	4
Connecting the Recirculation Cooling Unit.....	5
Connection to an External Exhaust System	5
2. Basic Operations	
Using the Controls.....	6
Front Panel Controls of Discotom-6	6
Front Panel Controls	7
Cooling Valve Positions	8
Cutting Position.....	8
Cleaning Position	8
Light Diode and Beep Signals.....	9
Fitting or Changing the Cut-off Wheel.....	10
Clamping the Workpiece	10
Starting/Stopping the Cutting Process.....	11
Automatic Cutting (AUTO mode).....	11
Manual Cutting	13
Combining Manual and Auto Operation	14
Cleaning	14
Flushing the Cutting Chamber.....	14

1. Getting Started

Checking the Contents

In the packing box you should find the following parts:

- 1 Fork spanner (24mm), for changing the cut-off wheel
- 1 Connector pipe for water outlet
- 1 Syringe with grease for Y-table spindle
- 1 Diagonal wrench for securing clamping devices
- 1 CD with Instruction Manuals

Recirculation Cooling Unit

- 1 Cooling water tank with lid
- 1 Recirculation pump
- 1 Recirculation sieve
- 1 Trolley
- 1 PVC hose ½", 2 m
- 1 Outlet hose, 2 m
- 2 Hose clamps
- 1 Large screw hose clamp
- 2 Disposable plastic inserts

Unpacking Discotom-6

Remove the screws underneath the pallet to free the Discotom-6.


Placing Discotom-6

Discotom-6 should be placed on a table that is able to carry a weight of min. 140 kg.


Discotom-6
Instruction Manual

**Getting Acquainted
with Discotom-6**

Take a moment to familiarise yourself with the location and names of the Discotom-6 components.

(Right side view)


- ① Cutting handle
- ② Protection guard, with window
- ③ Instrument panel
- ④ Emergency Stop
- ⑤ Flushing Gun
- ⑥ Cooling valve
- ⑦ Cutting table

- ⑧ Main power switch
- ⑨ Hours counter
- ⑩ Cut-off wheel locking handle

Discotom-6
Instruction Manual

Discotom-6 Rear View

- ① Upper back plate
- ② Electrical cable for power supply
- ③ Cut-off wheel locking handle
- ④ Reset button for recirculation pump
- ⑤ Power connection for recirculation pump
- ⑥ Exhaust hose
- ⑦ Water inlet
- ⑧ Water outlet
- ⑨ Hose for Flushing Gun

Noise Level

Approx. 67 dB (A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.

Power Supply

- First check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.
- Discotom-6 is factory mounted with an electric cable. Mount a plug on the cable:
Black and brown = phase
Yellow/green = earth

Direction of the Cut-off Wheel

Check that when the power is turned on, the cut-off wheel rotates in the direction indicated by the moulded arrow in the guard of Discotom-6. If the direction of rotation is incorrect, switch the two phases.

Connecting the Recirculation Cooling Unit

Connect the Recirculation cooling unit as follows:

- To connect the electric cable for the recirculation pump, remove the upper back plate on Discotom-6 and follow the instructions on the wiring diagram found on the inside of the back plate.
- Check that when the power is turned on, the pump rotates in the direction indicated by an arrow on the top of the pump. If not, re-check the connection of the electric cable.
- Replace the back plate.
- Clamp the inlet hose to the pump and attach the other end to the water inlet on Discotom-6 ⑦ .
- Attach the connector pipe to the water outlet ⑧, mount the large outlet hose and secure with the large hose clamp.
- Place the tank on the trolley.
- Place a disposable plastic insert in the tank and fold it over the rim.
- Fill the tank with 63 l water and 2 l Struers Additive. The water level should be 8-10 cm below the upper rim of the tank.
- Fit the tank lid and sieve.
- Roll the trolley into place and lead the outlet drain hose ⑨ into the cooling water tank. Adjust the lengths if necessary.

Reset Button

If the recirculation pump is overloaded or blocked, a thermal cut-out will switch the pump off. This safety cut-out must be reset by pushing the red button on the back of the Discotom-6 ④ before the machine can be restarted.

IMPORTANT

- ① Too high a level of coolant in the tank might damage the pump. To avoid this, place the disposable insert so that the pre-punched hole is in front of the overflow aperture in the tank.
- ② Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the Additive container). Remember to top up with Struers Additive each time you refill with water.


Connection to an External Exhaust System

Struers recommends the use of an exhaust as workpieces may emit harmful gasses or dust when cut. The exhaust system will also reduce the level of water condensation on the sides of the protection guard.

- On the back of Discotom-6 you will find a joint for a 50 mm dia. exhaust hose.
- Mount an exhaust hose from your local exhaust system onto the joint.

2. Basic Operations







Using the Controls Front Panel Controls of Discotom-6

Control name	Symbol
① Forwards button	▲
② Feed Speed Reduced indicator light.....	△
③ Backwards button	▼
④ Feed Speed potentiometer	
⑤ Load bar	
⑥ Emergency stop/overload indicator light	
⑦ Auto mode indicator light	Auto
⑧ Power indicator light	Power
⑨ Auto key	Auto
⑩ Start key	◇
⑪ Stop key	◊
⑫ Emergency stop button.....	

Discotom-6 Instruction Manual

Front Panel Controls


Name	Key	Function	Name	Key	Function
T MAIN SWITCH		The main power switch is located on the right hand side of Discotom-6. Turn clockwise to switch on the power.	① FORWARDS		Moves cutting table forward.
⑧ POWER INDICATOR		Lights when main power is switched on.	③ BACKWARDS		Moves cutting table backwards.
⑩ START		Starts the machine. The cut-off wheel starts rotating and the cooling water is turned on.	④ FEED SPEED		Regulates the feed speed of the cutting table. The speed can be set from 0.1-2.5 mm/sec.
⑪ STOP		Stops the machine. The cut-off wheel stops rotating and the cooling water is turned off.	② FEED SPEED REDUCED		Lights when the feed speed has been set too high and is reduced by the machine.
⑨ AUTO		Select or deselect automatic cutting mode	⑤ LOAD BAR		Reflects the present motor load status.
⑦ AUTO INDICATOR		Lights when automatic cutting mode selected.	⑫ EMERGENCY STOP		Push the red button to stop. Pull the red button to release.
			⑥ EMERG. STOP OVERLOAD		Lights when emergency stop has been pressed or the cutting motor is overloaded.

Cooling Valve Positions


The flushing gun, and a valve for cooling and cleaning, are located on the left hand side of the cutting machine.

During cutting the valve is turned clockwise (mark at the top) and the water is directed to cool the cut-off wheel. For cleaning the valve is turned counter-clockwise (mark at the bottom) to flush the cutting chamber.

Cutting Position

Cleaning Position

Discotom-6
Instruction Manual

Light Diode and Beep Signals

Discotom-6 is provided with various light diodes indicating the status of the machine. Beep signals will sound if a key is pressed when the function cannot be activated.

Status	Indication	Comments
Green POWER indicator on	The main power switch has been turned on.	
White AUTO indicator on	AUTO mode is active.	
🔔 Beep signal	You pressed AUTO, START  , FORWARDS  or BACKWARDS  while the motor was running.	
2 long 🔔 Beep signals	You pressed START  while the protection guard was open.	
3 long 🔔 Beep signals	You pressed START  for cutting while cooling valve isn't in cutting position.	Turn cooling valve to cutting position.
Red  indicator on	The emergency stop has been pressed.	Release the emergency stop to shut off the  indicator.
Red  indicator on and 🔔 beep signal	The emergency stop has been pressed and you tried to press another key.	
Yellow  FEED SPEED REDUCED indicator on	In AUTO mode: the feed speed is too high.	The machine automatically turns the feed speed down.
	An AUTO cut has finished. Feed speed reduced.	Turn the feed speed down or press any key.
Red  and yellow FEED SPEED REDUCED indicators on	Cutting in AUTO mode not possible as the cutting table does not move.	Press any key to stop the indicators. Material may be too hard for the cut-off wheel chosen or cutting table may be jammed. See Trouble Shooting section at back of manual.
Red  indicator and 1 red cell on LOAD BAR	Pump motor overloaded.	Stop the machine. Push the red button on the back of the machine to reset the safety cut-out.
Red  indicator and 2 red cells on LOAD BAR	Main motor overloaded. Process halted.	The indicators will disappear when the motor has cooled down.
Red  indicator, 2 red cells on LOAD BAR and 🔔 beep signal		You tried to press a key before the motor had cooled down.

Fitting or Changing the Cut-off Wheel

Note: The spindle for Discotom-6 is left-hand threaded. If Discotom-6 has been used for manual cutting, fasten the cutting handle with the black handle at the right hand side of the machine before changing the cut-off wheel.

- Press and hold down the black locking knob on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange and the old cut-off wheel.
- Mount the new cut-off wheel.

IMPORTANT

Conventional cut-off wheels based on $\text{Al}_2\text{O}_3/\text{SiC}$ abrasives should be placed between two cardboard washers, to protect the cut-off wheel. For maximum precision with diamond or CBN cut-off wheels, do not use cardboard washers.

- Mount the flange and the nut. Tighten carefully and release the locking knob.

Clamping the Workpiece

- Place the workpiece between the quick clamping devices and the back stops.
- Push the clamps towards the workpiece and lock the quick-clamping devices with the locking handle.

IMPORTANT

The workpiece must be firmly held in the quick-clamping devices.

- Carefully close the protection guard.

**Starting/Stopping
the Cutting Process**
*Automatic Cutting
(AUTO mode)*

When the machine is used in AUTO mode the cut-off wheel is stationary and the cutting table moves. In MANUAL mode the cutting table is stationary and the operator moves the cut-off wheel.

- Clamp the workpiece.
- Loosen the black handle situated on the right hand side of the cutting chamber.
- Lower the cut-off wheel by drawing the cutting handle downwards until the cut-off wheel is positioned ready to cut the workpiece. Fasten the black handle again.
- Press ▲ FORWARDS to position the cutting table and the workpiece just in front of the cut-off wheel.
- Close the protection guard.
- Set the desired feed speed
- Press the AUTO mode key to select AUTO mode.
- Press START ◇ to start the machine. The cut-off wheel starts rotating and the cooling water is turned on. The cutting table moves forward toward the cut-off wheel at the pre-set feed speed.
- If the machine has been stopped with the STOP ⊘ key, the cutting table will stay in position. To return the cutting table to its start position press BACKWARDS ▼ once. Struers recommends manually lifting the cut-off wheel away before the cutting table is moved back. Please note that the protection guard has to be closed in order to carry out this operation. If the protection guard has been opened press and hold the BACKWARDS key ▼ to move the cutting table to its start position.
- If the cutting table reaches the back of the cutting chamber during the cutting process, the cut-off wheel automatically stops and the cutting table returns to its start position.

■ To set another stop position:

- Press FORWARDS ▲ to move the cutting table to the desired position.
- Press and hold the STOP ⊖ key and then press once on FORWARDS ▲ to set this position as stop position. Then release the STOP ⊖ key.

When the cutting table now reaches this position during cutting, the cut-off wheel will stop and the cutting table will return to its start position.

The set stop position will be cancelled when the power is switched off.

Note: the cutting table can be positioned with FORWARDS ▲ and BACKWARDS ▼. When the keys are held down, the cutting table moves with increased speed of 20 mm / sec.

Reduced Feed Speed

If the feed speed is set too high, the yellow FEED SPEED REDUCED indicator (△) above the feed speed control flashes and the feed speed is automatically reduced by 10%. Discotom-6 can repeat this reduction up to 5 times.

Once the workpiece has been cut, the operator can reduce the feed speed using the feed speed control until the overload indicator stops flashing. Discotom-6 will then be able to cut another similar workpiece without feed speed reduction.

Manual Cutting

- Clamp the workpiece.
- Loosen the small black handle on the right side of the cutting chamber. The cutting handle is now easily moved up and down.
- Position the cutting table and the workpiece under the cut-off wheel with the FORWARDS ▲ and BACKWARDS ▼ keys
- Close the protection guard.
- Make sure that AUTO MODE is de-selected.
- Press START ◇ to start the machine. The cut-off wheel starts rotating and the cooling water is turned on.
- Pull down the cutting handle and let the cut-off wheel work itself into the workpiece. Increase the pressure and begin cutting.
- When the cut-off wheel has almost cut through the workpiece reduce the pressure.
- When the cut-off wheel has cut through the workpiece push the cutting handle back to its top position.
- Press the ⏻ STOP key to stop the machine

Remember


When cutting manually, the motor load bar should be used to monitor the force on the workpiece. Try to keep the motor load in the green or the beginning of the yellow area.

Combining Manual and Auto Operation

Manual and Automatic cutting mode may be used in combination. You may start in manual mode and cut into the workpiece. If you lock this position you can continue in automatic mode, and the workpiece will move towards the cut-off wheel.

Cleaning Flushing the Cutting Chamber

To ensure a longer lifetime for your Discotom-6, Struers strongly recommends that you clean the cutting chamber with the flushing gun every day.

- Press STOP ⏏ to stop the cut-off wheel and cooling water.
- Open the protection guard and release the workpiece from the quick-clamping devices.
- Take the flushing gun from the left hand side of the cutting machine and point it towards the bottom of the cutting chamber.
- Turn the cooling valve to cleaning position (see sketch at left)
- Clean the cutting chamber thoroughly.

Note:

Do not clean the protection guard with the flushing gun; use a damp cloth instead.

- Turn the cooling valve to cooling position.
- Leave the protection guard open to let the cutting chamber dry completely.

AVOID RUST!

Leave the protection guard open to let the cutting chamber dry completely.

Reference Guide

Table of Contents	Page
1. Advanced Operations	
Adjusting the Back Stop.....	16
Clamping Irregular Workpieces.....	16
Mounting an Adjustable Cutting Table.....	16
Optimising the Cutting Results.....	17
Explanation of Safety Factors.....	17
2. Accessories.....	18
3. Consumables.....	19
4. Trouble-Shooting	20
5. Maintenance	
Daily Service	23
Weekly Service.....	23
Checking the Recirculation Cooling Unit.....	23
Monthly service.....	23
Replacing the Cooling Water	23
Cleaning the Flushing Gun.....	23
Hours Counter	24
Maintenance of Cutting Tables.....	24
Lubricating the Spindle	24
Maintenance of Cut-off Wheels.....	25
Storing Bakelite Bonded Al ₂ O ₃ Cut-off Wheels	25
Maintenance of Diamond and CBN Cut-off Wheels	25
Maintenance of Clamping Devices.....	25
6. Technical Data	26

1. Advanced Operations

Adjusting the Back Stop

- Use a 13 mm spanner to loosen the back stop.
- Place the workpiece. Adjust the back stop so that the centre of the workpiece is under the centre of the cut-off wheel:
- Tighten the back stop.
- If you also need to adjust the clamping devices, loosen them and position them parallel to the back stop. Tighten the clamping devices.

Clamping Irregular Workpieces

Irregular workpieces without plane clamping surfaces must be clamped using special clamping tools, as the workpieces must not move during cutting. This could result in damage to the cut-off wheel or to the workpiece itself. For this purpose a vertical clamping device with one flat clamping shoe is available. (See Accessories).

Mounting an Adjustable Cutting Table

If you need to cut plane parallel sections, an adjustable table may be mounted on top of the existing cutting table. The adjustable table positions the workpiece at a right angle to the cut-off wheel and is movable up to 60mm. (See Accessories).

- Dismount the clamping tools, if mounted.
- Move the cutting handle to the top position and the cutting table to the back in order to create working space.
- Move the centre part of the adjustable cutting table to the left by turning the handle until the holes for the Allen screws are visible.
- Loosely mount two T-nuts underneath the adjustable cutting table using Allen screws.
- Slide the adjustable cutting table onto the cutting table:
 - Position the adjustable cutting table on the left-hand side of the cutting table with the spindle on the left.
 - Open the hinged protection guard to let the spindle pass.
 - The guide pins and T-nuts on the bottom of the adjustable cutting table are slid into the third T-slot from the left; this automatically aligns the adjustable cutting table.
- Tighten the Allen screws for the T-nuts.
- Mount clamping tools using the T-slots in the adjustable cutting table.

Optimising the Cutting Results

The following table shows answers to a number of common questions about how to achieve better quality cuts:

Question	Answer
How can I avoid discoloration or burning of the sample?	Reduce the feed speed
	Change the cut-off wheel, as the hardness of the present cut-off wheel may be inappropriate for the hardness of the workpiece*
How can I avoid burrs?	Use a softer cut-off wheel*
	Clamp the workpiece securely with the quick-clamping devices.
How can I prevent the cut-off wheel wearing out too quickly?	Reduce the feed speed or use a harder cut-off wheel*
How can I cut more quickly?	Place the workpiece in a position which allows the cut-off wheel to cut the smallest possible cross-section. Increase the feed speed.

*) Please refer to the Selection Guide in the Struers [Cut-off Wheels brochure](#).

Explanation of Safety Factors

The protection guard has a safety switch to prevent the cut-off wheel from starting while the guard is open. Furthermore, a locking mechanism prevents the operator from opening the protection guard until the cut-off wheel stops spinning.

2. Accessories

Specification	Code
Adjustable cutting table for cutting parallel sections	TRETA
Quick clamping device, left	TRELQ
Quick clamping device, right	TRERQ
Spring clamp, right	TRERI
Vertical clamping system with one flat clamping shoe	TREVS
Riser block for vertical clamping system (TREVS)	TREKS
Set of 4 swivel shoes for vertical clamping system (TREVS)	TREVI
Recirculation Cooling Unit, 65 l	TRECA
Replacement stainless steel bands 1 pc. 40 mm	TREXS
Replacement stainless steel bands For Adjustable Cutting Table 1 pc. 40 mm and 1 pc. 20mm	TRETT

3. Consumables

Cut-off Wheels

Please refer to the Selection Guide in the [Struers Cut-off Wheels brochure](#).

Other Consumables

Specification	Code
<i>Additive for Cooling Fluid</i> 1 l 5 l	ADDUN ADDFI
<i>Bacterial Cleaning Liquid</i> 0.5 l	BACTE
<i>Disposable Plastic Insert</i> for collecting waste material in the Recirculation cooling unit. Package with 10 pcs.	EXOSP

4. Trouble-Shooting

Error	Explanation	Action
Indications on the Front Panel		
FEED SPEED REDUCED $\triangle$ and ∇ indicator lights come on.	Emergency stop active.	Solve the problem causing the emergency stop, for example, by reducing the feed speed, or using a different cut-off wheel, then release the emergency stop.
	Cutting motor overloaded due to heavy use.	Open the protection guard and let the cutting motor cool down for about 5 - 10 min.
Machine Problems		
The cut-off wheel does not rotate.	The protection guard may be open.	Press the safety guard completely down. If this does not help, call a Struers service technician.
Chamber light does not work.	Lamp needs to be replaced.	Remove the plastic cap to gain access to the fluorescent lamp. Pull out the lamp and replace it.
Water leaking.	Leak in a Recirculation Cooling Unit hose.	Check the hose and tighten the hose clamp.
	Water overflow in the cooling water tank.	Remove the excess water in the tank.
Workpieces or cutting chamber rusty.	Insufficient additive in cooling fluid.	Use Struers Additive for cooling fluid in the cooling water, at the correct concentration. Check with a refractometer. Follow the instructions in the 'Maintenance' Section.
	The protection guard is left closed after use.	Leave the protection guard open to let the cutting chamber dry.
Protection guard window blurred.	Insufficient cleaning.	Clean with mild soapy water (do not use the flushing gun).
Quick-clamping device unable to hold the workpiece.	The quick-clamping device is not balanced.	Adjust the screw underneath the clamping column. Use a 3 mm Allen key.
	Clamp worn.	Call a Struers service technician.
Flushing gun blocked	Build up of dirt in the gun.	Clean flushing gun. Follow the instructions in the 'Maintenance' Section.

Discotom-6
Instruction Manual

Error	Explanation	Action
Cutting Problems		
Discoloration or burning of the workpiece.	The hardness of the cut-off wheel is inappropriate for the hardness / dimensions of the workpiece.	See Consumables section, under heading 'Cut-off Wheels'.
	Inadequate cooling.	- Check that there is sufficient water in the Recirculation cooling unit. - Check the cooling water hoses.
Unwanted burrs	Wheel is too hard.	See Consumables section, 'Cut-off Wheels'.
	Feed speed too high at the end of the operation.	Reduce the feed speed near the end of the operation.
	Lack of support for the workpiece.	If possible support the workpiece on both sides.
Cutting quality differs.	Cooling water hose clogged.	Clean the cooling water hose and the cooling tube. Check the water flow by turning the cooling valve to cleaning position.
Cut-off wheel breaks.	Incorrect mounting of the cut-off wheel.	- Check that the hole has the correct diameter. - Check cardboard washer on both sides of the cut-off wheels and replace if worn. The nut must be tightened properly.
	Incorrect clamping of the workpiece.	Make sure that only the left quick-clamping device is tight. The other device should only press lightly.
	Wheel is too hard.	See Consumables section, 'Cut-off Wheels'.
	Feed speed is set too high.	Reduce the feed speed.
	Inadequate cooling.	Check that there is enough water in the cooling unit. Check the cooling water hoses.

Discotom-6
Instruction Manual

Error	Explanation	Action
The cut-off wheel wears down too quickly.	The feed speed is too high.	Reduce the feed speed.
	Insufficient cooling.	- Check that there is enough water in the Recirculation cooling unit. - Check the cooling water hoses.
	The cut-off wheel is too soft for the task.	See Consumables section, 'Cut-off Wheels'.
The cut-off wheel does not cut through the workpiece.	Incorrect choice of cut-off wheel.	See Consumables section, 'Cut-off Wheels'.
	Cut-off wheel worn.	Replace the cut-off wheel.
	The cut-off wheel gets caught in the workpiece.	Support the workpiece and clamp it on both sides of the cut-off wheel so that the cut stays open.
The workpiece breaks when clamped.	The workpiece is brittle.	Place the workpiece between two polystyrene plates. Note! Always cut brittle workpieces very carefully.
The sample is corroded	The sample has been left in the cutting chamber for too long.	Remove the sample directly after cutting. Leave the protection guard open when you leave the machine.
	Insufficient additive for cooling fluid.	Use Struers Additive for cooling fluid in the cooling water at the correct concentration. Check with a refractometer. See Maintenance section.

5. Maintenance

Daily Service

- Clean the cutting chamber, especially the cutting table with the T-slots.
- Clean the protection guard window with a damp cloth. Do not use the flushing gun.
- Do **not** clean the lamp glass with alcohol. Use a damp cloth.
- Empty and clean the sieve in the Recirculation Cooling Unit.
- Refill the Recirculation Cooling Unit.

IMPORTANT

Always use the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to top up with Struers Additive each time you refill with water.

Weekly Service

Checking the Recirculation Cooling Unit

- Clean the cutting chamber thoroughly.

The cooling unit should be checked for cooling water after 8 hours use or at least every week. The unit must be refilled if the flushing pump cannot reach the cooling water. Remember to add Struers Additive for Cooling Fluid: One part of Additive for 33 parts of water.

To check the concentration of additive, use a refractometer. Concentration = Brix value. The concentration of additive should always lie between 2.7 and 3.3 %. Add Struers Additive for Cooling Fluid if the concentration is too low.

Monthly service

Replacing the Cooling Water

- Replace the cooling water in the Recirculation Cooling Unit at least once a month.
- Place a disposable plastic insert (ordered separately, please see section on *Consumables*) in the tank and fold it over the rim.
- **Please Note:** Too high a level of coolant in the tank might damage the pump. To avoid this, place the disposable insert so that the pre-punched hole is in front of the overflow aperture in the tank.
- Fill the tank with 63 l water and 2 l Struers Additive. The water level should be 8-10 cm below the upper rim of the tank.

Cleaning the Flushing Gun

- Remove flushing gun from the hose by squeezing the clamp.
- Soak in warm soapy water to loosen dirt and residue.
- Press the handle and clean the barrel with compressed air or a pipe cleaner. (The nozzle can be removed for better access).

Hours Counter

The hours counter at the right hand side of Discotom-6 records the number hours that the cutting motor has been in operation. To ensure trouble-free operation of Discotom-6, service provided by a Struers Service Technician is recommended for every 1500 hours of operation.

The hours counter can only be reset by a Struers Service Technician.

Maintenance of Cutting Tables *Lubricating the Spindle*

The stainless steel bands are also available as spare parts and should be replaced if they become worn or damaged.

To maintain the optimum performance of Discotom-6, lubricate the spindle of the cutting table at regular intervals (approx every 50 hours).

Check the service counter on the right hand side of the machine regularly to monitor the actual number of usage hours.

- Move the cutting table to its start position using the ▼ BACKWARDS key.
- Place the grease gun on the nipple in front of the spindle for the cutting table and push twice to grease the spindle.
- Note the date and number of service hours on the Maintenance Log table.

A syringe with grease for lubricating the Y-table spindle is supplied with the Discotom-6.

When all the grease has been used, refill the syringe with eg. Shell RETINAX Grease CS00 /Shell RETINAX Grease TL00.

Maintenance of Cut-off Wheels
Storing Bakelite Bonded Al_2O_3
Cut-off Wheels

These cut-off wheels are sensitive to humidity. Therefore, do not mix new, dry cut-off wheels with used humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.

Maintenance of Diamond and
CBN Cut-off Wheels

The precision of diamond and CBN cut-off wheels (and thus the cut) depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to a heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.
- Regular dressing of the cut-off wheel is also part of the general maintenance.

Maintenance of
Clamping Devices

Important

It is recommended to thoroughly clean and lubricate the Quick Clamping Devices at regular intervals.

6. Technical Data

Subject		Specifications	
		Metric/International	US
Cutting motor	Cutting power constant [S1] Cutting power intermittent [S3] Maximum power constant	2.2 - 3.0 kW 3.2 - 4.0 kW 4.4 - 6.0 kW	3.0 - 4.0 HP 4.3 - 5.4 HP 6.0 - 8.0 HP
	Rotational speed (running idle)	2850 rpm at 50Hz 3420 rpm at 60 Hz	2850 rpm at 50Hz 3420 rpm at 60 Hz
Cut-off Wheel	Diameter x thickness	ø250x1.5x32mm	9.8 x 0.06 x1.25"
	Max. cutting diameter: 80mm: Max. length of cut: 140mm (at a workpiece height of 65mm and width of 205mm) The max. cutting capacity might be less, depending on the workpiece material and the cut-off wheel	140 mm	5.6"
Cutting Table Dimensions	Width	386mm	15.19"
	Length	225mm	8.85"
Recirculation Cooling Unit	Approx. flow:	12 l/min	3.2 gallons/min
	Tank volume	65 l	17 gallons
Dimensions and Weight	Discotom-6		
	Width Depth Height Weight	794mm 640mm 458mm (123kg)	31.26" 25.20" 18.03" 272lbs
Environment	Safety standard	EU-directives 98/37/EEC - Safety of Machinery 73/23/EEC - Low Voltage Directive 89/336EEC, 92/31/EEC, 93/68/EF - EMC Directive Standards EN292-1, EN292-2, EN1050, EN60204-1 (IEC204-1), EN418, EN954-1, EN50081-1, EN50082-1, UL 508, NFPA70: 2002, FCC 47 Part 15 Class A	
	Dynamic Balancing	Max. admissible unbalance U_{per} according to ISO 1940/1, Balance Quality Grade G6.03: $U_{per} = 2 \text{ gmm/kg} = 110 \text{ gmm}$.	
Electrical Data	Noise level	Approx. 67 dB(A) during cutting, at a distance of 1.0 m/39.4" from the machine.	
	Voltage/frequency: 3 x 200-210V/50-60Hz 3 x 220-230V/50Hz 3 x 220-240V/60Hz 3 x 380-415V/50Hz 3 x 460-480V/60Hz	Max. Load: 21.3A 19.8A 19.8A 11.8A 11.8A	

Quick Reference Card

Clamping the Workpiece

- Place the workpiece between the clamps and the back stop.
- Push the clamps towards the workpiece and lock the quick-clamping device with the locking handle.

Auto Cutting

- Clamp the workpiece.
- Position the cut-off wheel so that it can cut the sample and lock it in place with the black cut-off wheel locking handle on the right of the machine.
- Position the cutting table with ▲ FORWARDS and ▼ BACKWARDS.
- Carefully close the protecting guard.
- Set feed speed and Press AUTO.
- Press START ◇ to begin cutting.
- After cutting, the cutting table automatically returns to its start position.
- Open the protection guard and release the clamping device.
- Remove the workpiece.

Manual Cutting

- Loosen the black handle on the right hand side of Discotom-6 so that the cut-off wheel moves freely
- Clamp the workpiece.
- Position the cutting table under the cut-off wheel.
- If AUTO mode is selected, de-select it.
- Press START ◇. Lower the cut-off wheel until it contacts the workpiece and apply pressure until the sample is cut.
- Keep the Load Bar in the green, or first part of the yellow area. Relax pressure near the end of the cut.
- Open the protection guard and release the clamping device.
- Remove the workpiece.

Changing the Cut-off Wheel

- Fasten the black handle to stop the cutting handle from moving.
- Press the black locking knob on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange and the old cut-off wheel.
- Mount the new cut-off wheel.
- Mount the flange and the nut. Tighten carefully.

Flushing the Cutting Chamber

- Pull out the flushing gun and point towards the bottom of the cutting chamber.
- Turn the cooling valve to flushing position.
- Clean the cutting chamber.
- Turn the cooling valve to cooling position.
- Always leave the protection guard open to let the cutting chamber dry.

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Date of Release 03.02.2006

Discotom-6
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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Spare Parts and Diagrams

Table of contents	Drawing
Discotom-6, complete	15670001F
Base and back plates, complete	15670015I
Base plate, complete.....	15670010G
Stepper motor and gear potentiometer	15670042C
Cutting spindle, complete.....	15670050E
Cutting wheel guard, complete.....	15110055A
Water valve, house	15110065D
Cover, complete	15670030E
Lamp, complete.....	15670012B
Cabinet, complete	15670020B
Front plate, complete	15670025B
Quick clamping device. Left	15110075C
Quick clamping device. Right.....	15110076A
Recirculation tank.....	14980068A
Recirculation pump	14980076D
	Diagrams
Wiring diagram (3 pages).....	15673450E
Circuit diagram	15673120D

Some of the drawings may contain position numbers not used in connection with this manual.

Discotom-6
Spare Parts and Diagrams

Spare part list for Discotom-6

Drawing	Pos.	Spare Part	Cat no:
15670001		Discotom-6, complete	
	210	O-ring 6.30-2.40 72 NBR 872	2IO24008
	335	50Hz coil	2MO01614
		60Hz coil	2MO01623
	465	W1 net cable 50Hz.	15113501
		W1 net cable 60Hz.	15113500
	470	Red Push Button For Thermal Relay.	2SA00301
15670015		Base and back plates, complete	
	0020	Motor 3x200V / 50Hz. 3x200-240V / 60 Hz.	2ME51140
		Motor 3x220-415V / 50 Hz	2ME51141
		Motor 3x460-480V / 60 Hz.	2ME51142
	0040	Key A 8X7X28	2ZF10828
	0190	Counter nut M25	2TD20026
	0200	Flex hose PMA PIST-17S.30	2NU31200
	0280	Cardboard flange Ø33/85	50000028
	0290	Flange	14150030
	0300	Washer for spindle	14150012
	0420	Edge list, black	2IP28459
	0430	Sealing band /31018 black	2IP30001
	0440	Sealing profile	2IP00203
	0560	Flushing gun	15670007
	0800	Thermo relay 0.18-0.27A	2FS10021
		Thermo relay 0.4-0.6A	2FS10023
		Thermo relay 0.6-0.9A	2FS10024
	0810	0.25 AT fuse, glass tube	2FU13500
15670010		Base plate, complete	
	20	Buffer AA-315 ø30x15 m. M8x10	2GS10315
	60	O-RING 50.00-4.00 NBR70	2IO40080
	390	Gas spring	15670046
	630	Transformer	2MT70711
	690	5.00A T Fuse glass tube	2FU14250
15670042		Stepper motor and gear potentiometer	
	0180	Stepper motor no.: 6600R174	2MI06674
	0260	Buffer ø30x40	2GS03040

Discotom-6
Spare Parts and Diagrams

Spare part list for Discotom-6

Drawing

15670050

Pos.	Spare Part	Cat no:
	Cutting spindle, complete	
10	Spindle shaft	15670110
20	Ball bearing, 6007 2RS1 ø35/ø62	2BK00110
30	Locking ring J62 DIN 472	2ZL20620
40	O-ring 31.42-1.78 NBR70	2IO17829
50	Seal ring	15110158
60+70	Seal ring	15112927
80	V-Ring V-40A	2IV10040
90	Lock nut	15110152
100	Cylinder pin, hardened. ST.6m6x14 D6325	2ZS03454
110	Spindle house painted	15679003
120	Spring	2GF51020
140	Lock ring	2ZL10350
150	Coupling	15670070
160	Disc M10	2ZA20010
170	ULS screw	2TR80812
180	Key 8x7x20	2ZF10820
185	Ring for sensor	15670411
186	Screw MSP M5x12	2TI10512
200	Namur MA516NF M5x16 sensor	2HQ00517

15110055

	Cutting wheel guard, complete	
20	Hose nipple 2601-12-3/8	2NF40088
30 → 90	Locking pin, complete	15112909
90	Grip M5	375MP085

15110065

	Water valve, house	
20	Inductive sensor	2HQ00023
30	O-ring 25.00-2.00	2IO20055

15670030

	Cover, complete	
20	Rubber curtain, sliced	15110665
30	Rubber curtain	15110666
60	Sealing profile EPMD (1M)	2IP21038

Discotom-6
Spare Parts and Diagrams

Spare part list for Discotom-6

Drawing

[15670012](#)

Pos.	Spare Part	Cat no:
	Lamp, complete	
0010	Lamp cover	15119001
0050	Bulb	2HG10827
0060	2uF 250V 50HZ condensator	2CR12301

[15670020](#)

	Cabinet, complete	
90	Seal profile 0,76 meter	2IP00205

[15670025](#)

	Front plate, complete	
20	Front plate foil	15670400
60	PCB Discotom-6 Contr. board	15673380

Accessories

Pipe for drain	15110997
Spanner for cutting wheel 24 mm	2GR00224
Socket spanner 13 mm	2GR00313
Manual CD	1067000
Grease gun (Shell Retinax CS00)	15670090
Grease for gun (Shell Retinax CS00)	2LS00000
Please specify amount	
Transportation crate	1511EM00

Wiring

A2 Brake circuit PCB	15670047
Contactor	2KM10640
Contactor K2-09A01 24V	2KM09024
Pot. Meter 5K 0,5W	2RP50120
W5 cable	15673505
W9 cable	15673509
w10 cable	15673510
W18 lamp cable	15673518

Spare part list for Discotom-6

Drawing


Pos.	Spare Part	Cat no:
	Various	
	Socket spanner (13 mm)	2GR00313
	Fork spanner (24 mm)	2GR00224


Vertical clamping system for 10 mm T-slots TREVS

15110075 Left / 15110076 Right


1	Operating key for TREVS	2GR27021
2	Spring for TREVS	2GR27023
3	Flat clamping shoe	2GR26150

Vertical clamping system for
10 mm T-slots (TREVS)











E	21-03-05	Pos 710 + 720 is change. + graphic update	MPE	21-03-05	CJE
D	10-01-05	Pos 290 new spring	MPE	10-01-05	CJE
C	03-12-2004	Pos. 240 fjernet	BMJ		
B	19-10-2004	Opdateret iht. stykl.	BMJ		
A	22-07-04		CJE		
Revision	Crea. date	Revision description	Draw. Init	Appr. date	Appr. Init
	04-07-04			04-07-04	
			Scale:	Format:	Tolerance: D5/ISO 2768-
			1:2.5	A1	Surface treat.:
ID: Description: 15670001 Discotom-6, complete					
Rev: E					





I	10-10-05	Pos. 910 is added	MPE	10-10-05	BMJ
H	09-08-2005	Pos. 115 added	BMJ		
G	04-07-05	Text for pos 20 is added the drawing	MPE	04-07-05	CJE
F	18-05-05	Pos 490 changed, graphic update	MPE	18-05-05	CJE
E	18-03-05	Pos 490, 500 & 510 changed	MPE	18-03-05	CJE
D	10-01-05	Pos 370 removed, Pos 480 & 490 changed	MPE	10-01-05	CJE
C	15-11-2004	Spulehoved erstattet af spulehane pos 560	MPE	15-11-04	CJE
B	18-10-2004	Pos. ner tilf. PF.15677401 tilf.	BMJ		
A					
Revision	Crea. date	Revision description	Draw. Init	Appr. date	Appr. Init
		Material:	Format:	Tolerance: DS/ISO 2768-	
		Scale: 1:2.5	A1	Surface treat.:	
ID:		Description:		Rev:	
15670015		Base and backplates assembled			


- 80 Bushings to be lubricated by Shell Helix 10W-40
- 130 Retinax CS00 (thin layer on the spindle)
- 270 to be filled with Shell Retinax Grease EP2
- 280 tightening torque = 35NM. Omnifit 50M
- 400 to be lubricated with Shell Retinax Grease EP2
- 740 Retinax Grease EP2 (on thread only)
- 790 To be seal with Omnifit 50H
- 80 Bßsninger smßres med Shell Helix 10W-40
- 130 Retinax CS00 (hele spindel smßres ind, tyndt lag)
- 270 fyldes med Shell Retinax Grease EP2
- 280 tilsp ndingsmoment = 35NM. Omnifit 50M
- 400 smßres med Shell Retinax Grease EP2
- 740 Shell Retinax Grease EP2 (kun p- gevind)
- 790 Tp tnes med Omnifit 50H

G	13-01-2006	Pos. 580 + 590(one time more) added.	BMJ	13-01-2006	JLI
F	25-10-05	Pos. 800 is added.	MPE	25-10-05	MPE
E	28-06-05	Pos. 790 is added.	MPE	28-06-05	CJE
D	18-05-05	Pos. 760 & 780 is added.	MPE	18-05-05	CJE
C	17-02-05	Pos 740+750 added. Pos 500/510/540-570 removed	JLI	17-02-05	JLI
B	26-10-04	Instructions for lubrication etc added	CJE	.	.
A	22-04-04		CJE	.	.
Revision	Crea. date	Revision description	Draw. Init	Appr. date	Appr. Init
		Material:	Scale: 1:5	Format: A1	Tolerance: DS/ISO 2768-Surface treat.:
		ID: 15670010 Base, assembled	Rev: G		

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


























Viewed at the back.

ALL WIRINGS 0.75mm², EXCEPT OTHERWISE MARKED.
 COLOR CODES: BK = Black, BN = Brown, RD = Red, OG = Orange, YE = Yellow, GN = Green, BU = Blue, VT = Violet, GY = Grey, WH = White, PK = Pink

WATER VALVE IN COOLING POSITION.


WATER VALVE IN SLICING POSITION.

M1 ROTATION DETECTION

RECIRCULATION IN OVERLOAD.

SAFETY GUARD OPEN.

EMERGENCY STOP ON

M1: CUTTING MOTOR

MANS SUPPLY:	Continuous (S1)	Non continuous (S1, 10% 12min)
3x200V / 50Hz	2.2 kW	3.2 kW
3x200 - 230V / 60Hz	4 HP	N/A
3x220 - 240V / 60Hz	4 HP	N/A
3x200 - 230V / 50Hz	2.2 kW	3.2 kW
3x380 - 415V / 50Hz	2.5 kW	3.2 kW
3x460 - 480V / 60Hz	4 HP	N/A

FT1: Thermal Overload Relay: SETTING VALUES

MANS SUPPLY:	FT1 SET =
3x200 - 230V / 50Hz	0.60 Amp
3x200 - 240V / 60Hz	0.53 Amp
3x380 - 415V / 50Hz	0.33 Amp
3x460 - 480V / 60Hz	0.30 Amp

S6 Bit Switch settings on Main Control Board

M1 Inom. current (from motor nameplate) *	Switch number	1	2	3	4
4.5 A	on	on	on	on	on
5.0 A	380-415V/50Hz	off	on	on	on
5.5 A	460-480V/50Hz	on	off	on	on
6.0 A		off	off	on	on
6.5 A		on	on	off	on
7.0 A		on	on	off	on
7.5 A		on	on	off	on
8.0 A	220-230V/50Hz	off	off	off	on
8.5 A		on	on	on	off
9.0 A		off	on	on	off
9.5 A		off	off	on	off
10.0 A		off	off	on	off
10.5 A		on	on	off	off
11.0 A	200-240V/50Hz	off	on	off	off
11.5 A	200V/50Hz	on	off	off	off
12.0 A		off	off	off	off

* Read motor nameplate carefully to decide if it shows S1- or S3-current values. Only if S1-current values is not available, use voltage variant column to decide correct switch settings.

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Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Denmark
Telephone +45 44 600 800

Herewith declares that

<i>Product Name:</i>	Discotom-6
<i>Type No:</i>	567
<i>Machine Type:</i>	Cut-off machine

is in conformity with the provisions of the following directives:

Safety of Machinery 98/37/EEC according to the following standard(s):
EN ISO 12100-1:2003, EN ISO 12100-2:2003/A1:1995, EN418:1992, EN1050:1996,
EN954-1:1996, EN60204-1:1997

EMC-Directive 89/336/EEC and 92/31/EEC according to the following standard(s):
EN50081-1:1992, EN50082-1:1997.

Low Voltage Directive 73/23/EEC according to the following standard(s):
EN60204-1:1997

Supplementary Information The equipment complies with the following standards:
NFPA70:2002, NFPA79:2002
FCC 47 CFR Part 15 Class A

The above has been declared according to the global method, module A

Date: 01.07.2006

Christian Skjold Heyde,
Vice President, R& D and Production, Struers A/S

Hersteller Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Danmark
Telefon +45 44 600 800

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<i>Produktname:</i>	Discotom-6
<i>Typennr.:</i>	567
<i>Maschinenart:</i>	Trennmaschine

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Christian Skjold Heyde,
Vice- President, R & D et Production, Struers A/S

Discotom-6

**Versatile automatic
and manual
cut-off machine in one**

- Automatic and manual cut-off machine in one
- Very high cutting capacity
- Constant feed speed principle ensures minimum deformation
- Automatic feed speed reduction for fast cutting of hard materials
- Quick positioning and freely selectable stop position functions
- Effective cooling system reduces the risk of burns
- Corrosion resistant cutting chamber
- Electric motor brake and safety lock ensure operator safety

Discotom-6 is a high performance table-top cut-off machine with a strong 3.7 kW (5.0 HP) motor. It is designed for cut-off wheels with a diameter of 250 mm (10").

Discotom-6 can be operated both automatically and manually, making it suitable for a wide range of applications.

Automatic Cutting made Easy

Automatic cutting mode is engaged by a simple press key. The cutting table with the clamped workpiece will move forward towards the cut-off wheel at the set feed speed. Before cutting, you can of course position the cut-off wheel to compensate for the wear of the cut-off wheel.

High Cutting Capacity

Manual mode for Fast and Uncomplicated Cutting

When operating the cut-off machine manually, the cutting table with the clamped workpiece is easily positioned under the cut-off wheel, and you pull the cutting handle to begin cutting.

Very High Cutting Capacity

You can cut workpieces with a diameter of up to 80 mm (3.1") on the Discotom-6. In automatic mode though, the cutting table is movable up to 200 mm (7.9") enabling cutting of very long samples. You can for instance cut a length of 140 mm (5.5") at a workpiece thickness of 65 mm (2.6").

The graph shows the cutting capacity in automatic cutting mode with a new cut-off wheel and flexible clamping tools. The real cutting capacity depends on the sample material

Quick positioning and Auto Stop

The cutting table can be moved with increased speed allowing quick positioning and set-up.


Additionally, a freely selectable stop position can be set where the cutting process will be stopped and the cutting table returns to its start position.

Automatic Feed Speed Reduction

When cutting automatically, you choose a feed speed from 0.1 to 2.5 mm/sec. If the feed speed is set too high and the machine cannot cut through the workpiece, Discotom-6 will automatically reduce the feed speed and continue cutting. This may repeat itself up to 5 times and prevents damage or burns on the sample surface.

Effective Cooling of the Cutting Area

A very important design detail is that the water cooling tubes are integrated into the cut-off wheel guard. The cooling water is applied to the cut-off wheel under high pressure and at a right angle. The cooling water thereby "breaks" the thin layer of air surrounding the spinning cut-off wheel, and the result is excellent cooling of the cutting area. Discotom-6 is available with a large 65 l recirculating cooling unit, ensuring that the cooling water does not overheat.

The adjustable cutting table enables parallel sectioning

Monitoring the Load on the Cutting Motor

The front panel includes a three-coloured LED-bar, indicating the load on the cutting motor.

In automatic cutting mode the LED-bar will show the load on the cutting motor and if the red area is reached the feed speed will be reduced. When cutting manually the LED-bar gives a good indication of how much force to put on the workpiece.

Corrosion resistant Cutting Chamber

The cutting chamber, made of corrosion resistant materials, is provided with a 225 mm x 386 mm (8.6" x 15.2") cutting table, which is easily moved up to 200 mm (7.9") by control keys. The cutting table is made of corrosion resistant cast aluminium and is fitted with a top surface

of stainless steel bands. These steel bands are easily exchanged if damaged or worn.

A flush hose makes daily cleaning quick and efficient, and a powerful lamp lights up the cutting chamber allowing you to monitor the cutting process through the large window in the protection guard. An exhaust system may be connected.

Quick and Handy Clamping of the Workpiece

The integrated cutting table has six 10 mm T-slots, facilitating the use of universal clamping tools.

Adding an Adjustable Cutting Table for Parallel Sectioning

An adjustable cutting table is available for mounting on the Discotom-6. It positions the workpiece at a right angle to the cut-off wheel, allowing parallel sections to be cut. The table can move up to 60 mm (2.4").

High Safety Level

Discotom-6 cannot be started while the protection guard is open. During cutting the protection guard is locked, eliminating any risk of injury. A built-in electric, dynamic motor brake stops the cutting motor within a few seconds, ensuring quick access to the cutting chamber after cutting.

For irregularly shaped workpieces the vertical clamping tools are employed

Quick-clamping devices enable easy and fast securing of workpieces

The hour counter registers the hours the cutting motor has been in operation

The flush hose facilitates cleaning

Technical Data

Cutting Motor	Cutting power constant [S1]	2.2 - 3.0 kW	3.0 - 4.0 HP
	Cutting power intermittent [S3, 15%]	3.2 - 4.0 kW	4.3 - 5.4 HP
	Maximum power	4.4 - 6.0 kW	6.0 - 8.0 HP
Cut-off Wheels	Cut-off wheel size: 250 mm x 1.5 mm x 32 mm (10" x 0.06" x 1.25") Rotational speed (running idle): 2850 rpm at 50 Hz, 3420 rpm at 60 Hz		
Cutting Table	Cutting table dimensions Width: 386 mm Length: 225 mm		
Dimensions and Weight	Discotom-6 without recirculating cooling unit Width: 794 mm (31.3") Depth: 640 mm (25.2") Height: 458 mm (18.1") Weight: 123 kg (272 lbs)		

Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact you local supplier for details)

Specifications

Discotom-6

Automatic cut-off machine. Automatic feed with electronic control of feed speed. Ability for manual operation. 3.7 kW (5.0 HP) motor. For 250 mm (10") cut-off wheels. Cutting table with 10 mm T-slots. Recirculation cooling unit (05766xxx), quick-clamping devices (05876909+05876910) or vertical clamping system (05116905) are ordered separately. Optional table unit (05116916) can be used.

Cat. no	Voltage
05676129	3 x 200 V / 50 Hz
05676130	3 x 200-210 V / 60 Hz CSA
05676135	3 x 220-230 V / 50 Hz
05676136	3 x 220-240 V / 60 Hz CSA
05676146	3 x 380-415 V / 50 Hz
05676147	3 x 380-415 V / 60 Hz
05676154	3 x 460-480 V / 60 Hz CSA

Accessories

Cooling System 3

With 50 l tank, small pump, Cooli-1, 50 l static filter with filter paper

Adjustable Cutting Table

For cutting plane parallel sections. It positions the workpiece at a right angle to the cut-off wheel. Movable up to 60 mm. With 10 mm T-slots and exchangeable stainless steel bands

Spring Clamp for 10 mm T-slots, Right

For retaining the cut-off specimen. Complete with back stop. To be mounted at the right hand side of the cut-off wheel

Vertical Clamping System for 10 mm T-slots

For clamping irregularly shaped workpieces on the cutting table. Clamping height adjustable up to 58 mm. Complete with operating key and one flat clamping shoe

Quick Clamping Tool for 10 mm T-slots, right

For securing the workpiece. To be mounted on the right side of the cut-off wheel. Complete with backstop. Max clamping width: 50 mm

Quick Clamping Tool for 10mm T-slots, left

For securing the workpiece. To be mounted on the left side of the cut-off wheel. Complete with backstop. Max clamping width: 50 mm

T-slot Turntable for 10mm T-slots

0-90 degree turntable for angular cutting. For use with various 10 mm T-slot clamping tools. Table height 42 mm.

Riser Block for 10 mm T-slots

For elevating Vertical Clamping System (05116905) when clamping high workpieces. To increase the clamping height by 60 mm

Swivel Shoes for Vertical Clamping System

Set of 4 multi-shaped swivel shoes for (05116905)

Replacement Stainless Steel Band

For cutting table on Discotom-6. 1 pc. of 40 mm x 225 mm

Cat. no.

056761xx

05766xxx

05116902

05116904

05116905

05876910

05876909

05876905

05116910

05116911

05676901

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

A wide range of cut-off wheels is available for Discotom-6. Please refer to the separate brochure.

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CANADA

Struers Ltd.
7275 West Credit Avenue
Mississauga, Ontario L5N 5M9
Phone +1 905-814-8855
Fax +1 905-814-1440
info@struers.com

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Struers Inc.